

## Application of Tape:

Prior to application check if the skin is dry and clean.

1. Apply tape smoothly.
  - Avoid creases or wrinkles.
2. Once tape has been applied, look to make sure surrounding skin's circulation is good.
  - Remove tape if too tight.

## Tape Removal:

Check for irritated or red skin.

1. Remove tape if any skin is irritated or red.
2. DO NOT pull or rip tape off aggressively.
3. Tape can be removed easily with water.



## CAUTION:

1. Avoid irritated, infected, or open wounds near skin to be taped.
2. Avoid taping skin if the patient has any history of hyper-sensitive skin or any pre-existing skin conditions.
3. Avoid any sensitive skin areas, or compromised circulation.



Manufactured for:

MedSource

8600 Shelby Court

Chanhassen, MN 55317 USA



Made in China

MSL-REG-QR-7168, REV. 0