ClearSafe® Safety IV Catheter

English

NOTE: U.S. Federal Law restricts this device to sale by or on the order of registered physician.
INTENDED USE: The MedSource ClearSafe® Safety IV Catheter is indicated to sample blood or administer fluids intravenously. The MedSource ClearSafe® Safety IV Catheter may be used for any patient population with consideration given to adequacy for the vascular anatomy appropriateness for the solution being administered and duration of the therapy.

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Product Description:

Each ClearSafe® Safety IV Catheter is comprised of the following components: protective needle cover, color-coded catheter hub, radiopaque catheter tube, teflon catheter holder, needle guide, medical grade stainless steel needle, needle hub, porous plug and gauge chamber.

The ClearSafe® Safety IV Catheter components combine for an ergonomic design that incorporates a safety feature. The safety feature of this device is intended to minimize accidental needle sticks. Please see the Instructions for Use section for further explanation of the safety feature.

The ClearSafe® Safety IV Catheter gauges are color coded for positive identification and range from 14G to 24G.

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Securely affix the catheter to the patient. Improper securement can result in loss of vascular access.

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Instructions for Use:

Due to the Risk of bloodborne pathogen exposure, follow Universal Precautions during placement, use, and removal of ClearSafe® Safety IV Catheter.

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	Product Ref. No.		Keep dry
	Batch Number		Fragile, handle with care
	Date Of Manufacturing		This way up
	Use By		Non-pyrogenic
	Do Not Reuse		Recyclable
	Sterilised by Ethylene Oxide Gas		Temperature limitation
	UDI		Importer
	Caution, Consult Accompanying Documents		Do not resterilize
	Consult Instructions for Use		Do Not Use if Package is Damaged
	Keep away from sunlight		Medical device
	Recyclable Packaging		Green Dot
	Rx only		Country of Manufacture

MedSource
8600 Shelby Court
Chanhassen MN, 55317 USA
Made in India

AW/IFU_AU, Rev. 03
Issue Date: 22.07.24
N. Code: HR/DEVICE/MFG/MD/2019/000135
MSL-REG-LQR-8604

ClearSafe Comfort® Safety IV Catheter

English

NOTE: U.S. Federal Law restricts this device to sale by or on the order of registered physician.
INTENDED USE: The MedSource ClearSafe Comfort® Safety IV Catheter is indicated to sample blood or administer fluids intravenously. The MedSource ClearSafe Comfort® Safety IV Catheter may be used for any patient population with consideration given to adequacy for the vascular anatomy appropriateness for the solution being administered and duration of the therapy.

The following instructions contain important information regarding the safe use of this product. Please read the Instruction for Use contents entirely before using this product. Failure to properly follow warnings and instructions could result in serious injury or death to the patient and/or clinician.

Product Description:

Each ClearSafe Comfort® Safety IV Catheter is comprised of the following components: protective needle cover, color-coded catheter hub, radiopaque catheter tube, teflon catheter holder, needle guide, medical grade stainless steel needle, needle hub, porous plug and gauge chamber.

The ClearSafe Comfort® Safety IV Catheter components combine for an ergonomic design that incorporates a safety feature. The safety feature of this device is intended to minimize accidental needle sticks. Please see the Instructions for Use section for further explanation of the safety feature.

The ClearSafe Comfort® Safety IV Catheter gauges are color coded for positive identification and range from 14G to 24G.

ORANGE = 14G GREY = 16G GREEN = 18G PINK = 20G BLUE = 22G YELLOW = 24G

Warnings:

Clinicians must follow Instructions for Use. In addition to following the Instructions for Use, clinicians must be trained to ensure proper usage of the device. Failure to do either of these could result in death or serious injury to the patient and/or clinician.

Best practice for avoiding an accidental needle stick is to use the device as described in the Instructions for Use and discard immediately into an approved sharps container.

Do not apply pressure during usage that causes the needle/catheter to flex.

In the unlikely event that the device does not function as described in the Instructions for Use, immediately dispose of the device in an approved sharps container.

Upon removal of the protective needle cover and prior to insertion, caution needs to be taken to ensure that the catheter does not come in contact with any foreign objects. Contact with a foreign object creates the potential for embolus. Never reinsert the introducer needle inside the catheter. Reinsertion threatens the integrity of the catheter and could lead to embolus.

Cautions:

This device in an unopened package is sterile, non-toxic and non-pyrogenic. Discard if open, wet, or damaged.

The materials used to manufacture this device do not contain natural rubber latex.

Not made with PVC. Not made with DEHP.

Radiopaque IV catheter.

This device is for single use only. Any attempt to re-sterilize and/or recondition the device will compromise the integrity of the device and will have an adverse affect on performance.

U.S. Federal Law restricts the sale of this device to or on the order of a registered physician.

Precautions:

This device is designed to reduce needle stick occurrences. Care must always be exercised when using this product. In addition to the Instructions for Use provided, it is strongly recommended that the clinician follow the guidelines put forth by the CDC and OSHA standards (USA) for bloodborne pathogens when using this product. During catheter insertion, maintenance, and removal, follow current "Guidelines for the Prevention of Intravascular Catheter-Related Infections", Centers for Disease Control and Prevention, Atlanta, GA (USA). Institutional policy and procedures and professional standard of care.

Securely affix the catheter to the patient. Improper securing can result in loss of vascular access.

A firm connection between the catheter hub and administration set is necessary to prevent leaking.

The introducer needle is longer than the catheter tube and will enter the vein or artery first. This can cause a blood flashback prior to the entrance of the catheter tube into the vein or artery. Advancing the catheter and needle together slightly will ensure intravascular access for the catheter. When advancing slightly, it is recommended to reduce the angle of approach in order to minimize the chance of puncturing the posterior wall of the vessel. In the event of an unsuccessful venipuncture, stabilize the catheter hub and actuate the device until it locks in the safety position. The safety position can be found by hearing a slight "click." Remove the catheter and properly dispose of both the spent catheter and locked Safety IV Catheter assembly.

Instructions for Use:

Due to the Risk of bloodborne pathogen exposure, follow Universal Precautions during placement, use, and removal of ClearSafe Comfort® Safety IV Catheter.

Clinicians must be trained in the practice of venipuncture and be aware of inherent dangers.

Keep unshielded needles away from body and fingers at all times.

Aseptic technique, proper skin preparation, and continued protection of the site are essential. Observe universal precautions on ALL patients.

- Remove needle cover in a straight outward motion and inspect catheter unit.
- Grasp device needle hub between thumb and fingers. Verify the release tab and needle bevel are in the up position.
- Insert needle into skin and vein and observe for blood flashback. The catheter is shorter than the needle. Blood flashback may occur before catheter tip is fully in the vein. It may be necessary to advance the catheter and needle slightly to achieve full entry into the vein. To avoid puncturing the posterior vein wall, lower the needle until it is parallel to the skin.
- Holding the flash chamber stationary, place finger on primary push off tab and thread desired catheter length into vein.
- Keeping the device stable, retract needle housing until it "clicks" into place. This is the indication that the needle is now locked inside the safety guard.
- Remove the needle guard by gently twisting and pulling it out of the catheter hub.
- Connect IV Administration line.
- Immediately discard unit into an approved sharps container.

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	Product Ref. No.		Keep dry
	Batch Number		Fragile, handle with care
	Date Of Manufacturing		This way up
	Use By		Non-pyrogenic
	Do Not Reuse		Recyclable
	Sterilised by Ethylene Oxide Gas		Temperature Limitation
	UDI		Importer
	Caution, Consult Accompanying Documents		Do not resterilize
	Consult Instructions for Use		Do Not Use if Package is Damage
	Keep away from sunlight		Medical device
	Recyclable Packaging		Green Dot
	Rx only		Country of Manufacture

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