


TrueSafe®

Push Button Style Safety IV Catheter



SPECIFICATIONS



- PTFE catheter material
- Push button style needle retraction
- Notched needle (sizes 20g - 24g)
- Ground stainless-steel needle
- Straight hub
- Disposable
- Not made with natural rubber latex
- Not made with DEHP
- Non-pyrogenic
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant
- Radiopaque




www.medsourcelabs.com

Item #	Description	UOM
MS-84414	TrueSafe® Safety IV Catheter, 14G x 1.75"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84416	TrueSafe® Safety IV Catheter, 16G x 1.77"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-844182	TrueSafe® Safety IV Catheter, 18G x 1.16"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84418	TrueSafe® Safety IV Catheter, 18G x 1.88"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-844202	TrueSafe® Safety IV Catheter, 20G x 1.16"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84420	TrueSafe® Safety IV Catheter, 20G x 1.88"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-844203	TrueSafe® Safety IV Catheter, 20G x 1.00"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84422	TrueSafe® Safety IV Catheter, 22G x 1.00"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84424	TrueSafe® Safety IV Catheter, 24G x 0.75"	50 EA/BX, 4 BX/CS, 200 EA/CS

TrueSafe Comfort®

Push Button Style Safety IV Catheter



SPECIFICATIONS



- VenaPUR® catheter material
- Push button style needle retraction
- Notched needle (sizes 20g - 24g)
- Ground stainless-steel needle
- Straight hub
- Disposable
- Not made with natural rubber latex
- Not made with DEHP
- Radiopaque
- Non-pyrogenic
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant




www.medsourcelabs.com


Item #	Description	UOM
MS-84314	TrueSafe Comfort®, 14G x 1.75"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84316	TrueSafe Comfort®, 16G x 1.77"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843182	TrueSafe Comfort®, 18G x 1.16"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84318	TrueSafe Comfort®, 18G x 1.88"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84320	TrueSafe Comfort®, 20G x 1.88"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843202	TrueSafe Comfort®, 20G x 1.16"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843203	TrueSafe Comfort®, 20G x 1.00"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84322	TrueSafe Comfort®, 22G x 1.00"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84324	TrueSafe Comfort®, 24G x 0.75"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843242	TrueSafe Comfort®, 24G x 0.56"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84326	TrueSafe Comfort®, 26G x 0.56"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843262	TrueSafe Comfort®, 26G x 0.75"	50 EA/BX, 4 BX/CS, 200 EA/CS

TrueSafe Comfort®

Push Button Style Safety IV Catheter



Audible Click Lock Safety Mechanism
Audible confirmation ensures the needle is safely encapsulated in the catheter body



SPECIFICATIONS



- VenaPUR® catheter material
- Push button style needle retraction
- Notched needle (sizes 20g - 24g)
- Ground stainless-steel needle
- Radiopaque
- Straight hub
- Disposable
- Not made with natural rubber latex
- Not made with DEHP
- Sterilized with ethylene oxide
- Non-pyrogenic
- ISO compliant color coding
- Needlestick Safety and Prevention Act (NSPA) of 2000 compliant

MEDSOURCELABS

www.medsourcelabs.com

Item #	Description	UOM
MS-84314	TrueSafe Comfort®, 14G x 1.75"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84316	TrueSafe Comfort®, 16G x 1.77"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843182	TrueSafe Comfort®, 18G x 1.16"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84318	TrueSafe Comfort®, 18G x 1.88"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84320	TrueSafe Comfort®, 20G x 1.88"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843202	TrueSafe Comfort®, 20G x 1.16"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843203	TrueSafe Comfort®, 20G x 1.00"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84322	TrueSafe Comfort®, 22G x 1.00"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84324	TrueSafe Comfort®, 24G x 0.75"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843242	TrueSafe Comfort®, 24G x 0.56"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-84326	TrueSafe Comfort®, 26G x 0.56"	50 EA/BX, 4 BX/CS, 200 EA/CS
MS-843262	TrueSafe Comfort®, 26G x 0.75"	50 EA/BX, 4 BX/CS, 200 EA/CS