



www.medsourcelabs.com

Reorder # MS-44000

QTY: 100

Adult Size

Blood Pressure Cuff Barrier

INSTRUCTIONS:

1. Slide barrier onto upper arm. Gather extra fabric at shoulder.
2. Apply Blood Pressure Cuff.
3. Fold extra fabric down to enclose Blood Pressure Cuff.
4. To remove, pull barrier and Blood Pressure Cuff as a unit.
Dispose the barrier after single use. Do not reuse.
Keep away from flame.



LOT:

CONTROL:



MedSource
8600 Shelby Court
Chanhassen, MN 55317
Made in China

